

Initial	comments	date	Rev
JW	Change distance of heatshrink tubing	7/10/19	Α

DRAWN BY			ALEPH OBJECTS® INCORPORATED		LULZBOT	
	West	2/27/19				
	CHECKED BY	DATE	TITLE			
			Y Bed Harness, Assembled			
	EL-HR0180		FILE NAME	SHEET	REV.	
			EL-HR0180_y_bed_harness_revA.odg	1	A	

AO part number	qty	UoM	Supplier	supplier pn
EL-WR0099	450	mm	Allcable	C0762A.41.10
EL-MS0247	20	mm	Electronics distributors corp	88-01870
PC-CN0003	1	pce	Digikey	WM2901-ND
EL-MS0059	5	pce	TTI global	16-02-0086
PC-CN0001	1	pce	Sager	50579402
PC-CN0119	1	pce	Digikey	WM2537-ND
EL-MS0205	5	pce	TTI global	16-02-0107
HD-TB0075	60	mm	McMaster	7856k76
EL-WR0157	430	mm	Wirecare	F6N0.25BK