

PETG DATA SHEET

This product has been GREENGUARD INDOOR AIR QUALITY CERTIFIED®. The GREENGUARD INDOOR AIR QUALITY CERTIFIED® Mark is a registered certification mark used under license through the GREENGUARD Environmental Institute (GEI). GEI is an industry independent, non-profit organization that oversees the GREENGUARD Certification Program. The GREENGUARD Certification Program is an industry independent, third-party testing program for low-emitting products and materials for indoor environments.

Mechanical Properties (Injection Molded), ATSM Method

Properties	Test Method	Typical Value
Specific Gravity	D 792	1.27
Water Absorption, 24h	D 570	0.1397
immersion	D 570	0.13%
Tensile Stress @ Break	D 638	28 MPa (4100 psi)
Tensile Stress @ Yield	D 638	50 MPa (7300 psi)
Elongation @ Break	D 638	130%
Tensile Modulus	D 638	2100 MPa (3.0 x 105 psi)
Flexural Modulus	D 790	2100 MPa (3.0 x 10⁵ psi)
Flexural Yield Strength	D 790	70 MPa (10200 psi)
Rockwell Hardness, R Scale	D 785	106
Izod Impact Strength, Notched		
@ 23°C (73°F)	D 256	101 J/m (1.9 ft·lbf/in.)
@ -40°C (-40°F)	D 256	37 J/m (0.7 ft lbf/in.)
Impact Strength, Unnotched		
@ 23°C (73°F)	D 4812	NB
@ -20°C (-4°F)	D 4812	NB
@ -30°C (-22°F)	D 4812	NB
@ -40°C (-40°F)	D 4812	NB
Impact Resistance (Puncture),		
Energy @ Max. Load	D 3763	
2.5-mm (0.100-in.) Thick		28 (21 ft (hf)
Plaques @ 23° C (73°F)	D 3763	20 5 (21 11 101)
	2 07 00	
2.5-mm (0.100-in.) Thick		41 J (30 ft·lbf)
Plaques, @ -40°C (-40°F)	D 3763	
3.2-mm (0.125-in.) Thick		33 J (24 ft ·lbf)
Plaques @ 23°C (73°F)	D 3763	

Mechanical Properties (Injection Molded), ISO Method

Properties	Test Method	Typical Value
Density	ISO 1183, Method D	1.27 g/cm3
Water Absorption, 24h	ISO 62	0.13%
Tensile Stress @ Break	ISO 527	28 MPa
Tensile Stress @ Yield	ISO 527	50 MPa
Elongation @ Break	ISO 527	100%
Tensile Modulus	ISO 527	2100 MPa
Flexural Modulus	ISO 178	2000 MPa
Flexural Yield Strength	ISO 178	68 MPa
Rockwell Hardness, R Scale	ISO 2039-2	109
Izod Impact Strength, Notched, Type 1 Specimen, Type A Notch		
@ 23°C	ISO 180	6.2 kJ/m2
@ -40°C	ISO 180	4.2 kJ/m2
Impact Strength, Unnotched, Type 1 Specimen ^a @ 23°C (73°F) @ -20°C (-4°F) @ -30°C (-22°F) @ -40°C (-40°F)	ISO 180 ISO 180 ISO 180 ISO 180	NB kJ/m2 NB kJ/m2 NB kJ/m2 NB kJ/m2
Impact Resistance (Puncture), Energy @ Max. Load ^b		
2.5-mm (0.100-in.) Thick Plaques, @ 23°C (73°F)	ISO 6603-2	40 J
2.5-mm (0.100-in.) Thick Plaques, @ -40°C (-40°F)	ISO 6603-2	35 J
3.2-mm (0.125-in.) Thick Plaques @ 23°C (73°F)	ISO 6603-2	44 J
3.2-mm (0.125-in.) Thick Plaques @ -40°C (40°F)	ISO 6603-2	36 J

Thermal Properties

Properties	Test Method	Typical Value
Deflection Temperature		
@ 0.455 MPa (66 psi)	D 648	70°C (158°F)
@ 1.82 MPa (264 psi)	D 648	64°C (147°F)
Vicat Softening Temperature	D 1525	85°C (185°F)
Thermal Conductivity	C 177	0.21 W/m·K (1.5 Btu in./h ft2.°F)
Glass Transition Temperature (Tg)	DSC	80°C (176°F)
Specific Heat		
@ 60°C (140°F)	DSC	1.30 kJ/kg⋅K (0.31 Btu/lb⋅°F)
@ 100°C (212°F)	DSC	1.76 kJ/kg·K (0.42 Btu/lb·°F)
@ 150°C (302°F)	DSC	1.88 kJ/kg⋅K (0.45 Btu/lb⋅°F)
@ 200°C (392°F)	DSC	1.97 kJ/kg⋅K (0.47 Btu/lb⋅°F)
@ 250°C (482°F)	DSC	2.05 kJ/kg ⋅K (0.49 Btu/lb ⋅°F)
Coefficient of Linear Thermal		5.1 x 10-5 /°C (mm/mm ·°C)
Expansion ^c	D 696	(2.8 x 10-5 /°F (in./in. ∘F))

Electrical Properties

Properties	Test Method	Typical Value
Dielectric Constant		
1kHz	D150	2.6
1MHz	D150	2.4
Dissipation Factor		
1kHz	D150	0.005
1MHz	D150	0.02
Arc Resistance	D 495	158 sec
Volume Resistivity	D 257	1015 ohm•cm
Surface Resistivity	D 257	1016 ohms/square
Dielectric Strength, Short Time, 500 V/sec rateof-rise	D 149	16 kV/mm (410 V/mil)

Note:

a Nonbreak as defined by ISO 180 with 4-mm specimens.

b Testing based on ISO 6603-2 using a striker diameter of 20 mm, a support and clamp diameter of 40 mm, and a velocity of 4.1 m/s.

c -30°C to 40°C (-22°F to 104°F)

Properties reported here are typical of average lots. Rigid.ink makes no representation that the material in any particular shipment will conform exactly to the values given.

Medical Disclaimer

The resin used in this filament has met selected FDAModified ISO-10993, Part 1 "Biological Evaluation of Medical Devices" tests with human tissue contact time of 30 days or less. The tests include: cytotoxicity, sensitization, irritation or intracutaneous reactivity, systemic toxicity (acute), subchronic toxicity (sub-acute), implantation, hemocompatibility. However it is the user who is responsible for the biological evaluation of the finished print for its intended application.

The suitability of a rigid.ink filament in a given end-use environment is dependent upon various conditions including, without limitation, chemical compatibility, temperature, part design, sterilization method, residual stresses, and external loads. It is the responsibility of the consumer to evaluate their final product under actual end-use requirements and to adequately advise and warn the end user where necessary.

www.rigid.ink

rigid.ink, Marston House, Walkers Court, Audby Lane, Wetherby, LS22 7FD

Reg. No. 9511774